


# INDIVIDUAL CANISTERS, FRAMES AND PLATFORM DIVIDERS

Canisters constructed of lightweight aluminum including locking handle and thumb hole. Frames constructed of corrosion resistant stainless steel including lifting handle. Custom sizes available upon request.

## INDIVIDUAL CANISTERS



### BAXTER FENWAL

<u>BAG NUMBER</u>	<u>PART NUMBER</u>	<u>DIMENSIONS A x B x C</u>
4R9951	ZC060H	3 <sup>5</sup> / <sub>8</sub> " x 6 <sup>1</sup> / <sub>4</sub> " x 3 <sup>3</sup> / <sub>8</sub> "
4R9953	ZC021	6" x 7 <sup>1</sup> / <sub>2</sub> " x 3 <sup>3</sup> / <sub>8</sub> "
4R9955	ZC022	6" x 9 <sup>3</sup> / <sub>8</sub> " x 3 <sup>3</sup> / <sub>8</sub> "
4R9957	ZC025	6" x 11" x 3 <sup>3</sup> / <sub>8</sub> "
4R9959	ZC023	6" x 14" x 3 <sup>3</sup> / <sub>8</sub> "
4R2111	ZC024	5 <sup>3</sup> / <sub>4</sub> " x 8" x 5 <sup>3</sup> / <sub>8</sub> "


### CHARTERMED

<u>BAG NUMBER</u>	<u>PART NUMBER</u>	<u>DIMENSIONS A x B x C</u>
7005-2	ZC004	7 <sup>1</sup> / <sub>4</sub> " x 9" x 3 <sup>3</sup> / <sub>8</sub> "
7039-2	ZC004	7 <sup>1</sup> / <sub>4</sub> " x 9" x 3 <sup>3</sup> / <sub>8</sub> "
7303-2	ZC004	7 <sup>1</sup> / <sub>4</sub> " x 9" x 3 <sup>3</sup> / <sub>8</sub> "
7400-2	ZC005	9" x 9" x 3 <sup>3</sup> / <sub>8</sub> "
7410-2	ZC005	9" x 9" x 3 <sup>3</sup> / <sub>8</sub> "
7403-2	ZC005	9" x 9" x 3 <sup>3</sup> / <sub>8</sub> "
SPECIAL	ZC002	11 <sup>1</sup> / <sub>2</sub> " x 11" x 3 <sup>3</sup> / <sub>8</sub> "
PEDIATRIC	ZC001	5 <sup>3</sup> / <sub>4</sub> " x 9 <sup>1</sup> / <sub>8</sub> " x 7 <sup>1</sup> / <sub>16</sub> "

### GAMBRO

<u>BAG NUMBER</u>	<u>PART NUMBER</u>	<u>DIMENSIONS A x B x C</u>
DF-200	ZC011	8" x 6 <sup>3</sup> / <sub>16</sub> " x 1 <sup>3</sup> / <sub>16</sub> "
DF-700	ZC012	12" x 6 <sup>3</sup> / <sub>16</sub> " x 1 <sup>3</sup> / <sub>16</sub> "
DF-1000	ZC013	14 <sup>1</sup> / <sub>4</sub> " x 6 <sup>3</sup> / <sub>16</sub> " x 1 <sup>3</sup> / <sub>16</sub> "
DF-1200	ZC014	16 <sup>1</sup> / <sub>2</sub> " x 6 <sup>3</sup> / <sub>16</sub> " x 1 <sup>3</sup> / <sub>16</sub> "


## INDIVIDUAL FRAMES



<u>FRAME</u>	<u>PART NUMBER</u>
2 PLACE	ZF2
3 PLACE	ZF3
4 PLACE	ZF4
5 PLACE	ZF5
6 PLACE	ZF6
7 PLACE	ZF7
8 PLACE	ZF8

NOTE: To order, add canister part number after frame part number.  
Example: ZF4-ZC021

## INDIVIDUAL PLATFORM DIVIDERS



<u>FREEZER DIAMETER</u>	<u>PART NUMBER</u>
18 INCH	ZD18
19 INCH	ZD19
20 INCH	ZD20
31 INCH	ZD31
39 INCH	ZD39

NOTE: To order, add canister part number after platform divider part number. Example: ZD31-ZC021

# COMPLETE SYSTEMS CANISTERS, FRAMES AND PLATFORM DIVIDERS

Systems consist of complete set of canisters, frames plus platform divider. Canisters constructed of lightweight aluminum including locking handle and thumb hole. Frames including lifting handle and platform dividers are constructed of corrosion resistant stainless steel. Custom sizes available upon request.

## 20 INCH FREEZER

### BAXTER FENWAL

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
4R9951	434	62 (7 PL)	ZS200	27.75"
4R9953	224	56 (4 PL)	ZS201	26"
4R9955	164	41 (4 PL)	ZS202	26"
4R9957	136	34 (4 PL)	ZS203	26"
4R9959	88	22 (4 PL)	ZS204	26"

### CHARTERMED

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
7005-2	129	43 (3 PL)	ZS205	24"
7039-2	129	43 (3 PL)	ZS205	24"
7303-2	129	43 (3 PL)	ZS205	24"
7400-2	129	43 (3 PL)	ZS206	27.75"
7410-2	129	43 (3 PL)	ZS206	27.75"
7403-2	129	43 (3 PL)	ZS206	27.75"
SPECIAL	70	35 (2 PL)	ZS207	25"

### GAMBRO

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
DF-200	108	27 (4 PL)	ZS208	27"
DF-700	64	16 (4 PL)	ZS209	27"

## 39 INCH FREEZER

### BAXTER FENWAL

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
4R9951	1932	276 (7 PL)	ZS900	27.75"
4R9953	912	228 (4 PL)	ZS901	26"
4R9955	752	188 (4 PL)	ZS902	26"
4R9957	608	152 (4 PL)	ZS903	26"
4R9959	456	114 (4 PL)	ZS904	26"

### CHARTERMED

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
7005-2	552	184 (3 PL)	ZS905	24"
7039-2	552	184 (3 PL)	ZS905	24"
7303-2	552	184 (3 PL)	ZS905	24"
7400-2	552	184 (3 PL)	ZS906	27.75"
7410-2	552	184 (3 PL)	ZS906	27.75"
7403-2	552	184 (3 PL)	ZS906	27.75"
SPECIAL	304	152 (2 PL)	ZS907	25"

### GAMBRO

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
DF-200	464	116 (4 PL)	ZS908	27"
DF-700	320	80 (4 PL)	ZS909	27"

## 31 INCH FREEZER

### BAXTER FENWAL


<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
4R9951	1120	160 (7 PL)	ZS100	27.75"
4R9953	560	140 (4 PL)	ZS101	26"
4R9955	440	110 (4 PL)	ZS102	26"
4R9957	352	88 (4 PL)	ZS103	26"
4R9959	256	64 (4 PL)	ZS104	26"

### CHARTERMED

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
7005-2	336	112 (3 PL)	ZS105	24"
7039-2	336	112 (3 PL)	ZS105	24"
7303-2	336	112 (3 PL)	ZS105	24"
7400-2	336	112 (3 PL)	ZS106	27.75"
7410-2	336	112 (3 PL)	ZS106	27.75"
7403-2	336	112 (3 PL)	ZS106	27.75"
SPECIAL	180	90 (3 PL)	ZS107	27"

### GAMBRO

<u>BAG NUMBER</u>	<u>CANISTER QUANTITY</u>	<u>FRAME QUANTITY</u>	<u>PART NUMBER</u>	<u>OAH</u>
DF-200	288	72 (4 PL)	ZS108	27"
DF-700	192	48 (4 PL)	ZS109	27"





## Viewing and Printing Tips for Acrobat PDF Documents

Note: the interface and specific techniques may vary depending on the version of Acrobat Reader that you are using. Adobe Acrobat Reader is a free program. You may download the latest version from Adobe.com. To see complete instructions for Acrobat, use the HELP button at the top of Acrobat.

### BETTER VIEWING

The clarity of a PDF document on your screen may vary depending on the size you choose to view the document at.

You can read dimensions better and generally view a document better by zooming in. To zoom in, click the Zoom In button on the top toolbar (the + sign). Or, choose Zoom In from the View menu. You can decrease the size of the document with the Zoom Out tool (the - sign)

To zoom in on a specific area, click the Zoom button on the toolbar, then click in the area on the document which you want to enlarge. Zooming out works the same way.

You can also view a document at a specific zoom level. Select Zoom To from the View menu. The zoom dialog box appears. Type the zoom percentage you desire in the text box and click OK. The document will be displayed at your specified magnification.

### PRINTING

Regardless of the magnification and clarity of the document you are viewing, printing the document will generally result in a very clear image of text and product dimensions. These Acrobat documents are optimized to print clearly on your local computer printer.

To print the document to your local computer printer, either choose Print from the File menu, click the Print button on the File toolbar, or press Ctrl P (Command P on the Mac). The Print dialog box appears.

You may select the number of copies you want to print in the Print dialog box.

If you don't want to print the entire document, enter the page range (the starting and ending page numbers of the pages to print) in the appropriate text boxes of the Print dialog box.

Click OK (Windows) or Print (Mac)

For additional printing instructions, use the Help button on the top toolbar of Acrobat Reader.

Note: if you need to print this catalog in part or in its entirety through a commercial printer, contact Custom BioGenic Systems for an Adobe Acrobat PDF version suitable for high-end commercial printing.