

Ready-to-use cleaner  
for system and capillaries.

CASY®clean



▶ CASY®clean

*Performance starts here*

Schärfe System sets the standard:  
both in technology and supplies.



## CASY<sup>®</sup>clean: ready-to-use cleaner for system and capillaries

CASY supplies are specially matched to the standardized CASY devices. Only the exclusive and regular use of CASY supplies guarantees long-term standardization of measurement results, as well as the long-term stability and optimal functionality of the equipment.



CASYclean is a ready-to-use cleaner specified by Schärfe System for use with the CASY Cell Counter + Analyser Systems.

### Description

CASYclean is a foam-free, alkaline, odorless and colorless cleaning solution for the system and measurement capillaries. CASYclean is non toxic and non flammable.

### Properties

Its regular use protects device and measurement capillaries against biological deposits. Through its emulsifying effect, CASYclean prevents the internal liquid containing system and the measurement capillaries from blockages and encrustations.

### Packaging

CASYclean is available from Schärfe System direct in a pack of 3 x 500 ml bottles.

### Shelf-life and storage

The best-before date is indicated on the product packaging. CASYclean should be stored in a cool, dry place, protected from light.

### Waste disposal

Disposal must be effected in accordance with the cell culture laboratory regulations.

**SCHÄRFE**   
**SYSTEM**

Schärfe System GmbH  
Krämerstraße 22  
D-72764 Reutlingen  
Phone: +49 (0) 7121 / 387 86-0  
Fax: +49 (0) 7121 / 387 86-99  
mail@CASY-Technology.com  
www.CASY-Technology.com

 Put us to the test.