












## ScanSpeed Centrifuge Product Guide

Centrifuges		Low-Speed Centrifuges		Multi-Purpose Centrifuges				High-Speed	Microcentrifuges			
												
Refrigeration					•		•		•	•		
Max, Force												
Fixed Angle Rotors	Max, 2,000 x g		•									
	Max, 2,700 x g	•										
	Max, 12,300 x g											•
	Max, 16,500 x g			•	•							
	Max, 21,200 x g					•					•	
	Max, 25,900 x g						•					
	Max, 27,200 x g								•	•		
	Max, 32,300 x g								•			
	Max, 54,111 x g							•				
Swing Rotors	Max, 2,800 x g	•										
	Max, 3,100 x g						•	•				
	Max, 3,600 x g					•	•					
	Max, 5,300 x g			•	•							
Sample Containers												
Microtubes	PCR Strips(x8)									•	•	•
	0,2/0,5 ml			•	•	•	•	•	•	•	•	•
	1,5/2,0 ml			•	•	•	•	•	•	•	•	•
Tubes	3/4/5/10/15 ml	•	•	•	•	•	•	•	•			
	50 ml	•		•	•	•	•	•	•			
	100 ml			•	•	•	•	•	•			
	250 ml			•	•	•	•	•	•			
	500 ml					•	•	•	•			
	750 ml					•	•					
	1,000 ml							•				
Microplates	PCR plate	•		•	•	•	•	•	•			
	Micro titer	•		•	•	•	•	•	•			
	Cell culture	•		•	•	•	•	•	•			
	Deep well			•	•	•	•	•	•			

**a/s ninolab**

DIN LEVERANDØR AF BASISUDSTYR TIL LABORATORIER

Tel 44911007 • [www.ninolab.dk](http://www.ninolab.dk)